

Product Datasheet

SLC6A8 peptide (orb217939)

Description	The peptide is used to block anti-SLC6A8 antibody reactivity.
Reactivity	Human, Mouse, Rabbit, Rat
Tested Applications	Blocking
Target	SLC6A8
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214583.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q8VBW1, P48029, P28570
Entrez	50690, 6535, 102857
Expiration Date	12 months from date of receipt.