

Product Datasheet

SLC8A1 peptide (orb217940)

Description	The peptide is used to block anti-NCX1 antibody reactivity.
Reactivity	Human, Mouse, Rabbit, Rat
Tested Applications	Blocking
Target	SLC8A1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214584.
Purity	> 85%
Source	Synthetic
Uniprot ID	P70414, Q01728, P32418
Entrez	29715, 6546, 20541
Expiration Date	12 months from date of receipt.