

Product Datasheet

SLC16A2 peptide (orb217942)

Description	The peptide is used to block anti-MCT8 antibody reactivity.
Reactivity	Human, Mouse, Rabbit, Rat
Tested Applications	Blocking
Target	SLC16A2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214586.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q8K1P8, P36021, O70324
Entrez	20502, 6567, 259248
Expiration Date	12 months from date of receipt.