

Product Datasheet

SLC22A2 peptide (orb217945)

Description	The peptide is used to block anti-OCT2 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	SLC22A2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214589.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9R0W2 , O15244 , O70577
Entrez	20518 , 6582 , 29503
Expiration Date	12 months from date of receipt.