

Product Datasheet

SMARCA2 peptide (orb217946)

Description	SMARCA2 peptide
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	SMARCA2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214590.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q6DIC0, P51531
Entrez	6595, 67155
Expiration Date	12 months from date of receipt.