

Product Datasheet

SMPD1 peptide (orb217947)

Description	The peptide is used to block anti-Acid Sphingomyelinase antibody reactivity.
Reactivity	Human, Mouse, Rat, Sheep
Tested Applications	Blocking
Target	SMPD1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214591.
Purity	> 85%
Source	Synthetic
Uniprot ID	P17405, Q04519
Entrez	6609, 20597
Expiration Date	12 months from date of receipt.