

Product Datasheet

SNAP25 peptide (orb217950)

Description	The peptide is used to block anti-SNAP25 antibody reactivity.
Reactivity	Human, Mouse, Porcine, Primate, Rat, Zebrafish
Tested Applications	Blocking
Target	SNAP25
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214594.
Purity	> 85%
Source	Synthetic
Uniprot ID	P60880, P60881, P60879
Entrez	6616, 25012, 20614
Expiration Date	12 months from date of receipt.