

Product Datasheet

alpha Synuclein peptide (orb217952)

Description	alpha Synuclein peptide
Reactivity	Human, Mouse, Porcine, Rat, Sheep
Tested Applications	Blocking
Target	SNCA
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214596.
Purity	> 85%
Source	Synthetic
Uniprot ID	P37840, P37377, O55042
Entrez	29219, 6622, 20617
Expiration Date	12 months from date of receipt.