

Product Datasheet

SNCG peptide (orb217953)

Description	The peptide is used to block anti-Gamma-synuclein antibody reactivity.
Reactivity	Human, Primate
Tested Applications	Blocking
Target	SNCG
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214597.
Purity	> 85%
Source	Synthetic
Uniprot ID	076070
Entrez	6623
Expiration Date	12 months from date of receipt.