

Product Datasheet

SOX9 peptide (orb217957)

Description	The peptide is used to block anti-SOX9
Reactivity	Human, Mouse, Porcine, Rat
Tested Applications	Blocking
Target	SOX9
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214601.
Purity	> 85%
Source	Synthetic
Uniprot ID	P48436, Q04887
Entrez	6662, 20682
Expiration Date	12 months from date of receipt.