

Product Datasheet

SOX11 peptide (orb217958)

Description	The peptide is used to block anti-SOX11 antibody reactivity.
Reactivity	Human, Mouse, Rat, Zebrafish
Tested Applications	Blocking
Target	SOX11
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P0C1G9, Q7M6Y2, P35716
Entrez	6664, 20666, 84046
Expiration Date	12 months from date of receipt.