

Product Datasheet

SRC peptide (orb217962)

Description	The peptide is used to block anti-c-SRC
Reactivity	Human, Mouse, Rat, Zebrafish
Tested Applications	Blocking
Target	SRC
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214606.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9WUD9, P05480, P12931
Entrez	6714, 20779, 83805
Expiration Date	12 months from date of receipt.