

Product Datasheet

SRF peptide (orb217965)

Description	The peptide is used to block anti-SRF antibody reactivity.
Reactivity	Human, Mouse
Tested Applications	Blocking
Target	SRF
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214609.
Purity	> 85%
Source	Synthetic
Uniprot ID	P11831, Q9JM73
Entrez	20807, 6722
Expiration Date	12 months from date of receipt.