

Product Datasheet

STAT3 peptide (orb217976)

Description	The peptide is used to block anti-STAT3 antibody reactivity.
Reactivity	Human, Mouse, Porcine, Rat
Tested Applications	Blocking
Target	STAT3
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214620.
Purity	> 85%
Source	Synthetic
Uniprot ID	P42227, P40763, P52631
Entrez	25125, 6774, 20848
Expiration Date	12 months from date of receipt.