

Product Datasheet

STAT5A peptide (orb217980)

Description	The peptide is used to block anti-STAT5A antibody reactivity.
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	STAT5A
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214624.
Purity	> 85%
Source	Synthetic
Uniprot ID	P42229, P42230
Entrez	6776, 20850
Expiration Date	12 months from date of receipt.