

Product Datasheet

STAT6 peptide (orb217982)

Description	STAT6 peptide
Reactivity	Human, Mouse, Porcine, Primate, Rabbit
Tested Applications	Blocking
Target	STAT6
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214626.
Purity	> 85%
Source	Synthetic
Uniprot ID	P52633, P42226
Entrez	6778, 20852
Expiration Date	12 months from date of receipt.