

Product Datasheet

STRN peptide (orb217986)

Description	The peptide is used to block anti-Striatin antibody reactivity.
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	STRN
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214630.
Purity	> 85%
Source	Synthetic
Uniprot ID	O55106, P70483, O43815
Entrez	6801, 268980, 29149
Expiration Date	12 months from date of receipt.