

Product Datasheet

SULT2A1 peptide (orb217987)

Description	The peptide is used to block anti-DHEA-ST antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	SULT2A1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214631.
Purity	> 85%
Source	Synthetic
Uniprot ID	P15709, Q06520, P52843
Entrez	6822, 20859, 24912
Expiration Date	12 months from date of receipt.