

Product Datasheet

SYK peptide (orb217991)

Description	The peptide is used to block anti-SYK
Reactivity	Human, Mouse, Porcine, Rabbit, Rat
Tested Applications	Blocking
Target	SYK
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214635.
Purity	> 85%
Source	Synthetic
Uniprot ID	P48025, Q64725, P43405
Entrez	20963, 6850, 25155
Expiration Date	12 months from date of receipt.