

Product Datasheet

SYP peptide (orb217995)

Description	The peptide is used to block anti-Synaptophysin antibody reactivity.
Reactivity	Human, Mouse, Porcine, Primate, Rat
Tested Applications	Blocking
Target	SYP
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214639.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q62277, P08247, P07825
Entrez	24804, 6855, 20977
Expiration Date	12 months from date of receipt.