

## Product Datasheet

### ADAM17 peptide (orb217996)

|                            |                                                                          |
|----------------------------|--------------------------------------------------------------------------|
| <b>Description</b>         | The peptide is used to block anti-CD156b antibody reactivity.            |
| <b>Reactivity</b>          | Human, Mouse, Rat                                                        |
| <b>Tested Applications</b> | Blocking                                                                 |
| <b>Target</b>              | ADAM17                                                                   |
| <b>Form/Appearance</b>     | Lyophilized powder                                                       |
| <b>Storage</b>             | Shipped at 4°C. Store at -20°C for one year.                             |
| <b>Note</b>                | For research use only                                                    |
| <b>Application notes</b>   | This peptide was used as immunogen to produce orb214640.                 |
| <b>Purity</b>              | > 85%                                                                    |
| <b>Source</b>              | Synthetic                                                                |
| <b>Uniprot ID</b>          | <a href="#">Q9Z0F8</a> , <a href="#">P78536</a> , <a href="#">Q9Z1K9</a> |
| <b>Entrez</b>              | <a href="#">6868</a> , <a href="#">57027</a> , <a href="#">11491</a>     |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                          |