

Product Datasheet

TADA2A peptide (orb217997)

Description	The peptide is used to block anti-ADA2A antibody reactivity.
Reactivity	Human, Mouse, Porcine, Primate, Rabbit, Rat
Tested Applications	Blocking
Target	TADA2A
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214641.
Purity	> 85%
Source	Synthetic
Uniprot ID	O75478, Q8CHV6, Q6AYE3
Entrez	360581, 217031, 6871
Expiration Date	12 months from date of receipt.