

Product Datasheet

TFAP2A peptide (orb218009)

Description	The peptide is used to block anti-AP2 alpha/beta antibody reactivity.
Reactivity	Human, Mouse, Porcine, Rat, Sheep
Tested Applications	Blocking
Target	TFAP2A; TFAP2B
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214653.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q92481, Q61313, P34056, P58197, P05549
Entrez	21418, 7020, 21419
Expiration Date	12 months from date of receipt.