

Product Datasheet

TFF2 peptide (orb218014)

Description	The peptide is used to block anti-Spasmolysin antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	TFF2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214658.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q03404, Q03403, Q09030
Entrez	7032, 116592
Expiration Date	12 months from date of receipt.