

Product Datasheet

TGF Beta 1 peptide (orb218018)

Description	TGF Beta 1 peptide
Reactivity	Human, Mouse, Porcine, Rat
Tested Applications	Blocking
Target	TGFB2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214662.
Purity	> 85%
Source	Synthetic
Uniprot ID	P61812, P27090, Q07257
Entrez	81809, 21808, 7042
Expiration Date	12 months from date of receipt.