

Product Datasheet

TGFB3 peptide (orb218019)

Description	The peptide is used to block anti-TGF beta 3 antibody reactivity.
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	TGFB3
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214663.
Purity	> 85%
Source	Synthetic
Uniprot ID	P17125, Q07258, P10600
Entrez	7043, 25717
Expiration Date	12 months from date of receipt.