

Product Datasheet

THRA peptide (orb218026)

Description	THRA peptide
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	THRB
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214670.
Purity	> 85%
Source	Synthetic
Uniprot ID	P10828
Entrez	7068
Expiration Date	12 months from date of receipt.