

Product Datasheet

TIMP3 peptide (orb218029)

Description	The peptide is used to block anti-TIMP3 antibody reactivity.
Reactivity	Human, Mouse, Porcine, Rabbit, Rat
Tested Applications	Blocking
Target	TIMP3
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214673.
Purity	> 85%
Source	Synthetic
Uniprot ID	P48032, P35625, P39876
Entrez	21859, 7078
Expiration Date	12 months from date of receipt.