

Product Datasheet

ICAM5 peptide (orb218032)

Description	The peptide is used to block anti-ICAM-5 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	ICAM5
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214676.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q60625, Q9UMF0
Entrez	15898, 7087
Expiration Date	12 months from date of receipt.