

Product Datasheet

TNF alpha peptide (orb218037)

Description	TNF alpha peptide
Reactivity	Human, Mouse, Porcine, Primate, Rat
Tested Applications	Blocking
Target	TNF
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214681.
Purity	> 85%
Source	Synthetic
Uniprot ID	P06804, P16599, P01375
Entrez	7124, 24835, 21926
Expiration Date	12 months from date of receipt.