

Product Datasheet

TNFR1 peptide (orb218040)

Description	TNFR1 peptide
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	TNFRSF1A
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214684.
Purity	> 85%
Source	Synthetic
Uniprot ID	P25118, P22934, P19438
Entrez	7132, 21937, 25625
Expiration Date	12 months from date of receipt.