

Product Datasheet

PGP9.5 peptide (orb218071)

Description	The peptide is used to block anti-PGP9.5 antibody reactivity.
Reactivity	Human, Mouse, Porcine, Primate, Rat
Tested Applications	Blocking
Target	UCHL1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214715.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q00981 , Q9R0P9 , P09936
Entrez	22223 , 29545 , 7345
Expiration Date	12 months from date of receipt.