

Product Datasheet

VASP peptide (orb218077)

Description	The peptide is used to block anti-VASP antibody reactivity.
Reactivity	Human, Primate, Rat
Tested Applications	Blocking
Target	VASP
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214721.
Purity	> 85%
Source	Synthetic
Uniprot ID	P50552
Entrez	7408
Expiration Date	12 months from date of receipt.