

Product Datasheet

VCP peptide (orb218081)

Description	The peptide is used to block anti-VCP antibody reactivity.
Reactivity	Human, Mouse, Porcine, Rat, Zebrafish
Tested Applications	Blocking
Target	VCP
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214725.
Purity	> 85%
Source	Synthetic
Uniprot ID	P46462, P55072, Q01853
Entrez	7415, 269523, 116643
Expiration Date	12 months from date of receipt.