

Product Datasheet

VHL peptide (orb218085)

Description	The peptide is used to block anti-Von Hippel Lindau antibody reactivity.
Reactivity	Human, Rat
Tested Applications	Blocking
Target	VHL
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214729.
Purity	> 85%
Source	Synthetic
Uniprot ID	P40337, Q64259
Entrez	7428, 24874
Expiration Date	12 months from date of receipt.