

Product Datasheet

EZR peptide (orb218087)

Description	The peptide is used to block anti-Ezrin antibody reactivity.
Reactivity	Human, Mouse, Porcine, Rabbit, Rat
Tested Applications	Blocking
Target	EZR
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214731.
Purity	> 85%
Source	Synthetic
Uniprot ID	P26040, P31977, P15311
Entrez	7430, 54319, 22350
Expiration Date	12 months from date of receipt.