

Product Datasheet

YWHAG peptide (orb218105)

Description	The peptide is used to block anti-14-3-3 gamma antibody reactivity.
Reactivity	Human, Mouse, Rat, Zebrafish
Tested Applications	Blocking
Target	YWHAG
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214749.
Purity	> 85%
Source	Synthetic
Uniprot ID	P61982, P61983, P61981
Entrez	22628, 7532, 56010
Expiration Date	12 months from date of receipt.