

Product Datasheet

14-3-3 peptide (orb218107)

Description	The peptide is used to block anti-14-3-3 zeta antibody reactivity.
Reactivity	Human, Mouse, Porcine, Rat, Sheep, Zebrafish
Tested Applications	Blocking
Target	YWHAZ
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214751.
Purity	> 85%
Source	Synthetic
Uniprot ID	P63104, P63101, P63102
Entrez	7534, 22631, 25578
Expiration Date	12 months from date of receipt.