

Product Datasheet

ZAP70 peptide (orb218109)

Description	The peptide is used to block anti-ZAP70
Reactivity	Human, Mouse, Porcine, Primate, Rabbit, Rat
Tested Applications	Blocking
Target	ZAP70
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214753.
Purity	> 85%
Source	Synthetic
Uniprot ID	P43404, P43403
Entrez	7535, 22637
Expiration Date	12 months from date of receipt.