

Product Datasheet

CNBP peptide (orb218110)

Description	The peptide is used to block anti-ZNF9 antibody reactivity.
Reactivity	Human, Mouse, Porcine, Rat
Tested Applications	Blocking
Target	CNBP
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214754.
Purity	> 85%
Source	Synthetic
Uniprot ID	P53996, P62634, P62633
Entrez	12785, 64530, 7555
Expiration Date	12 months from date of receipt.