

Product Datasheet

CXCR4 peptide (orb218117)

Description	CXCR4 peptide
Reactivity	Human, Mouse, Porcine, Primate, Rat
Tested Applications	Blocking
Target	CXCR4
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214761.
Purity	> 85%
Source	Synthetic
Uniprot ID	O08565, P61073, P70658
Entrez	12767, 7852
Expiration Date	12 months from date of receipt.