

Product Datasheet

SLC39A7 peptide (orb218119)

Description	SLC39A7 peptide
Reactivity	Human, Mouse, Porcine, Rat
Tested Applications	Blocking
Target	SLC39A7
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214763.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q92504, Q31125
Entrez	7922, 14977
Expiration Date	12 months from date of receipt.