

Product Datasheet

MTERF peptide (orb218121)

Description	The peptide is used to block anti-mTERF antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	MTERF
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214765.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q8CHZ9, Q99551, Q9EPI8
Entrez	545725, 7978, 85261
Expiration Date	12 months from date of receipt.