

Product Datasheet

FZD6 peptide (orb218132)

Description	The peptide is used to block anti-Frizzled 6 antibody reactivity.
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	FZD6
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214776.
Purity	> 85%
Source	Synthetic
Uniprot ID	O60353, Q61089
Entrez	14368, 8323
Expiration Date	12 months from date of receipt.