

Product Datasheet

HIST1H3A peptide (orb218137)

Description	The peptide is used to block anti-Histone H3
Reactivity	Human, Mouse, Porcine, Rat, Zebrafish
Tested Applications	Blocking
Target	HIST1H3A
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214781.
Purity	> 85%
Source	Synthetic
Uniprot ID	P68431, P84244, P84243, P84245, P68433, Q71DI3, Q6LED0, P84228
Entrez	15078, 8358, 8352, 8351, 8355, 8357, 8968, 291159, 100361558, 8350, 8353, 8354, 15077, 8356, 319152
Expiration Date	12 months from date of receipt.