

Product Datasheet

SLC25A11 peptide (orb218138)

Description	The peptide is used to block anti-SLC25A11 antibody reactivity.
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	SLC25A11
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214782.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q02978 , Q9CR62 , P97700
Entrez	67863 , 8402 , 64201
Expiration Date	12 months from date of receipt.