

Product Datasheet

CUL2 peptide (orb218141)

Description	The peptide is used to block anti-Cullin 2 antibody reactivity.
Reactivity	Human, Mouse, Porcine, Rat
Tested Applications	Blocking
Target	CUL2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214785.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9D4H8, Q13617
Entrez	8453, 71745
Expiration Date	12 months from date of receipt.