

Product Datasheet

RGS5 peptide (orb218146)

Description	The peptide is used to block anti-RGS5 antibody reactivity.
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	RGS5
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214790.
Purity	> 85%
Source	Synthetic
Uniprot ID	P49800, O15539, O08850
Entrez	8490, 54294, 19737
Expiration Date	12 months from date of receipt.