

Product Datasheet

SLC43A1 peptide (orb218148)

Description	The peptide is used to block anti-SLC43A1 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	SLC43A1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214792.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q8BSM7, O75387
Entrez	8501, 72401
Expiration Date	12 months from date of receipt.