

Product Datasheet

PEX3 peptide (orb218149)

Description	The peptide is used to block anti-PEX3 antibody reactivity.
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	PEX3
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214793.
Purity	> 85%
Source	Synthetic
Uniprot ID	P56589 , Q9JJK4 , Q9QXY9
Entrez	83519 , 56535 , 8504
Expiration Date	12 months from date of receipt.